Channel classification instruction v3

Cheng-Wei Shih, Chia-Ming Kuo

National Central University, Taiwan



- Create a folder, make sure it has
 - run.sh
 - summary_plot.c
 - calibration_ana_code_multi.c
 - check_chip_prototypeMaximam_new.c
 - All ".dat" files (not .root file), 10 runs is enough.

Example

```
[[5202011@chip01 INTT_multi_run_v2_test3]$ ls

NCU_fphx_raw_module_268_20210610-1132_0.dat

NCU_fphx_raw_module_268_20210610-1137_0.dat

NCU_fphx_raw_module_268_20210610-1137_0.dat

NCU_fphx_raw_module_268_20210610-1143_0.dat

NCU_fphx_raw_module_268_20210610-1148_0.dat

NCU_fphx_raw_module_268_20210610-1148_0.dat

NCU_fphx_raw_module_268_20210610-125_0.dat

NCU_fphx_raw_module_268_20210610-125_0.dat

NCU_fphx_raw_module_268_20210610-125_0.dat

NCU_fphx_raw_module_268_20210610-123_0.dat

NCU_fphx_raw_module_268_20210610-1231_0.dat

NCU_fphx_raw_module_268_20210610-1204_0.dat

NCU_fphx_raw_module_268_20210610-1241_0.dat

NCU_fphx_raw_module_268_20210610-1241_0.dat

NCU_fphx_raw_module_268_20210610-1241_0.dat
```

Root version: 6.01

Before run the code, please modify "run.sh"

```
older_directio<mark>n≧"</mark>/home/5202011/INTT
                                                                                  Folder direction
 number_of_file=1
merge_file_name="aaa_test_summary" Name of final merged root file
 nodule_ID=2
rm multi_run_status.txt
rm $merge_file_name.root
echo 1
                                         Module ID
sleep 15
ls *.dat > dat_file.txt
sleep 15
let number_for_final=number_of_file-1
for seed in $(seq 0 $number_for_final)
       cp check_chip_prototypeMaximam_new.c check_chip_prototypeMaximam_new_copy.c
       sed -i "s/data_index/${seed}/g" check_chip_prototypeMaximam_new_copy.c
       root -l -b -g check_chip_prototypeMaximam_new_copy.c\($module_ID\)
       rm check_chip_prototypeMaximam_new_copy.c
       sleep 15
ls *.root > total_file.txt
sleep 15
for seed in $(seq 0 $number_for_final)
   cp calibration_ana_code_multi.c calibration_ana_code_multi_copy.c
   sed -i "s/data_index/${seed}/g" calibration_ana_code_multi_copy.c
   root -l -b -g calibration_ana_code_multi_copy.c\(\"$folder_direction\",$module_ID,true,false,0,true,false,false,true\)
   rm calibration_ana_code_multi_copy.c
                                                      Variables descriptions are in next slide
   sleep 15
root -l -b -q summary_plot.c\(\shumber_of_file,\"\sfolder_direction\",\"\sherge_file_name\"\)
```

Variable 3 & 6: if "false" -> no plots created, can be faster

```
#Variable of calibration_ana_code_multi.c :
 ariable 1 : TSting, folder direction
 ′ariable 2 : int,
                   port_ID,
                   output the adc-ampl plot for each channel (should be true)
 ′ariable 3 : bool,
                   original unbond channel check
                                                      (should be false)
                   overall ampl noise level check
                                                  (0 can be good)
                   output offset ampl distribution plot for each channel (should be true)
                   cout unbonded channel status @ without bias run. (should be false)
                   cout wider gaus width channel @ with bias run. (should be false)
#Variable 8 : bool.
                                                           (should be true)
                   output multi run status.txt
```

- After modification, please run "run.sh" → ". run.sh"
 - It takes ~ 15 mins to finish 20 root files.


- Two files and a lot folders will be created after the run
 - Two files: ppb2_l2_summary.root & multi_run_status.txt
 - One .root file will be created in each folder of each run to record the status of the problem channels.
- 3 final summary plots and un-functional channels status will be created and printed, examples are shown in next slide.
 - 3 final summary plots:
 - Summary_noise.pdf
 - Summary_entry.pdf
 - Summary_average_entry.pdf

Final overall results

Channels failed in noise criteria

```
Noise channel found, chip: 2
                              channel : 0
                                             failed times : 13/20
                                                                    ratio : 0.65
Noise channel found, chip: 7
                              channel: 12
                                             failed times : 20/20
                                                                    ratio: 1
 Noise channel found, chip : 19 channel : 36
                                             failed times : 20/20
                                                                    ratio : 1
 Noise channel, chip: 2 channel: 0
 failed in file index : 2
                              file name: fphx_raw_20210428-1320_0.root, gaus width: 46.9425
 failed in file index : 3
                              file name: fphx_raw_20210428-1321_0.root, gaus width: 78.376
 failed in file index : 4
                              file name: fphx_raw_20210428-1322_0.root, gaus width: 4.47558
                              file name : fphx_raw_20210428-1324_0.root, gaus width : 4.59855
 failed in file index : 5
 failed in file index : 6
                              file name: fphx_raw_20210428-1331_0.root, gaus width: 75.7756
 failed in file index: 7
                              file name: fphx_raw_20210428-1332_0.root, gaus width: 32.8639
 failed in file index : 8
                              file name: fphx_raw_20210428-1333_0.root, gaus width: 38.3274
 failed in file index : 11
                              file name: fphx_raw_20210428-1654_0.root, gaus width: 147.086
 failed in file index : 12
                              file name: fphx_raw_20210428-1656_0.root, gaus width: 127.588
 failed in file index : 14
                              file name: fphx_raw_20210428-1659_0.root, gaus width: 4.8683
                              file name: fphx_raw_20210428-1700_0.root, gaus width: 4.15262
 failed in file index : 15
 failed in file index: 17
                              file name: fphx_raw_20210428-1703_0.root, gaus width: 4.41939
 failed in file index: 19
                              file name: fphx_raw_20210428-1706_0.root, gaus width: 33.9221
Channels failed in entry criteria
Weird entries found, chip: 2 channel: 0
                                             failed times : 12/20
                                                                    ratio: 0.6
Weird entries found, chip: 19 channel: 113 failed times: 20/20
                                                                   ratio : 1
 Bad entries channel, chip: 2 channel: 0
 bad in file index : 2 file name : fphx_raw_20210428-1320_0.root, entries : 4410
 bad in file index: 3 file name: fphx_raw_20210428-1321_0.root, entries: 5820
 bad in file index : 5 file name : fphx_raw_20210428-1324_0.root, entries : 645
 bad in file index: 6 file name: fphx_raw_20210428-1331_0.root, entries: 7557
 bad in file index: 7 file name: fphx_raw_20210428-1332_0.root, entries: 2104
 bad in file index: 8 file name: fphx_raw_20210428-1333_0.root, entries: 4062
 bad in file index: 11 file name: fphx_raw_20210428-1654_0.root, entries: 20615
 bad in file index: 12 file name: fphx_raw_20210428-1656_0.root, entries: 17784
 bad in file index: 14 file name: fphx_raw_20210428-1659_0.root, entries: 569
 bad in file index: 15 file name: fphx_raw_20210428-1700_0.root, entries: 422
 bad in file index: 17 file name: fphx_raw_20210428-1703_0.root, entries: 467
 bad in file index: 19 file name: fphx_raw_20210428-1706_0.root, entries: 2125
             ----- Final summary -----
# of noise channels : 3
# of bad entry chan : 2
# of double counting: 1
Total bad channels: 4
```

Principle: the less entries in the plot, the better performance it is



Summary_average_entry.pdf

Criteria of bad channel

- # of input files > 10 files.
- Single channel is tested by 2 criteria :
 - Entry:
 - If channel with entry < 280 OR entry > 400.
 - If it happens more than 2 times in all calibration runs -> Bad.
 - Noise:
 - If channel with fit gaussian width > 4.
 - If it happens more than 2 times in all calibration runs -> Bad.
- If channel fails in both criteria, only counts 1.
 - No double counting.

Plot descriptions

There are a lot of plots in each folder created by ". run.sh". Here I introduce some plots I frequently check


- chipX_detail_ampladc.pdf
 - Ampl ADC response for single channel
 - 128 pages for 128 channels




- Ampl width distribution after offset, single channel.
- 128 pages for 128 channels




- Overall gaus width status of half-ladder.
- Entries of plot : 128*26 = 3328



X axis : ampl, Y axis : ADC






Plot descriptions

There are a lot of plots in each folder created by ". run.sh". Here I introduce some plots I frequently check


- ampl_adc_width_detial.pdf
 - Overall gaus width status of each chip
 - 26 pages in total
- channel_entries.pdf
 - Distribution of # of events of each channel
 - Entries: 26*128 = 3328



of data points: 128


X axis: Channel,


Y axis: Gaussian width



Attention!

- When you check "chipX_detail_amploffset.pdf", if the distribution is not reasonable (For example : double or triple peaks, left plot). Please let me know.
- It is a bug, and has been solved already.







Bug case, The amtpl-adc distribution is good, but it has a bug when filling the event in TH1F

Algorithm introduction

- For each channel of each chip :
 - Center line: mean of events in "Block 0" (ADC==0)
 - Center of the rest blocks are panned to center line.
 - Amount of movement: Mean of each block center line
 - Each event is filled in TH1F, and fit with gaussian.




X axis : amp, Y axis : ADC




Algorithm introduction

Update of my algorithm, 2 criteria: noise and entries

Gaus width > 4 will be shown in right plot



Algorithm introduction



Back up